

Pulsed Interleaved Excitation (PIE)

PIE stands for **P**ulsed **I**nterleaved **E**xcitation: a technique of synchronising several pulsed lasers (or other suitable excitation sources). The laser pulses are separated on a ns time scale to allow simultaneous recording of the time response of a sample molecule due to each individual laser.

Copyright of this document belongs to PicoQuant GmbH. No parts of it may be reproduced, translated or transferred to third parties without written permission of PicoQuant GmbH. All information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuraciesor omissions. Specifi cations and external appearances are subject to change without notice.



PicoQuant GmbH Rudower Chaussee 29 (IGZ) 12489 Berlin Germany P +49-(0)30-1208820-89 F +49-(0)30-1208820-90 info@picoquant.com www.picoquant.com

© PicoQuant GmbH 1/1