

Software

Commercial

Microscopy

[SymPhoTime64](#)

A complete software suite for acquisition and analysis of time resolved fluorescence data. SymPhoTime features a powerful scripting language that allows customisation of existing analysis and even implementation of novel analysis techniques.

[Phasor Analysis](#) Phasor Analysis is currently not included in [SymPhoTime64](#). However, you can use [Globals](#) developed by the [Laboratory for Fluorescence Dynamics](#) to analyse FLIM data via the Phasor approach.

Spectroscopy

[EasyTau](#)

EasyTau is the systems software of PicoQuants time resolved spectrometer [FluoTime 300](#). EasyTau features a scripting language to automate more complex experiments. Examples of EasyTau scripts can be found on github: <https://github.com/PicoQuant/EasyTau-Scripts>

[FluoFit](#)

Open Source

Please note that these are external projects for which PicoQuant takes no responsibility and provides no support. In case of problems please contact the respective authors. Also note that such projects may not be maintained continuously. In particular, please check if they use current versions of file formats and/or programming libraries before you invest in larger development work based on such projects.

If you develop improved or new code, please consider contributing your code back to the community.

The compilation below is probably not complete. If you are the maintainer of a software that deals with TCSPC data please edit this page and add your software listing. Alternatively you can also use <https://support.picoquant.com> to contact PicoQuant support to submit your software project.

Data Analysis

- [Demo Code for reading PicoQuants PTU File format](#)
- [A Python library to read PTU files](#)
- [ptuT3toT2](#) - a program that converts ptu files from picoquant from T3 mode to T2 mode
- [Extensible Time-tag Analyzer](#)

- [exciton diffusion \(Matlab\)](#)
- [Reading PicoQuants PTU File formats to MatLab - 100x faster than the Demo Code](#)
- [Researchgate group Open Computational Tools for Single-Molecule Spectroscopy](#)
- PAM is a software package for quantitative analysis of fluorescence microscopy and spectroscopy data, with a focus on experiments using pulsed interleaved excitation. <http://www.cup.uni-muenchen.de/pc/lamb/software/pam.html>
- ImageJ Import Filter for pt3 files: http://imagejdocu.tudor.lu/doku.php?id=plugin:inputoutput:picoquant_.pt3_image_reader:start
- ImageJ import of ptu/pt3 image files (+intensity/average lifetime stacks): https://github.com/ekatruxha/PTU_Reader
- **DecayFit - Fluorescence Decay Analysis** <http://www.fluortools.com/software/decayfit>
- **AniFit - Anisotropy Decay Analysis** <http://www.fluortools.com/software/anifit>
- [PyCorrFit](#)
- [PyCorrelate](#)
- QuickFit 3.0: <http://www.dkfz.de/Macromol/quickfit/>
- [FCS viewer](#)
- FCS Bulk Correlation Software accepting pt3 and ptu files
[FCS Point Correlator](#)
- [FLIMFit](#)
- **Photon-Tools** Python utilities for working with photon timestamp data from fluorescence spectroscopy experiments <https://github.com/bgamari/photon-tools>
- **Photon Correlation** A suite of command-line tools for performing calculations on photon arrival data, as produced by hardware such as the PicoQuant PicoHarp and HydraHarp: https://github.com/tsbischof/photon_correlation
- **libpicoquant:** Tools for reading PicoQuant file formats: <https://github.com/tsbischof/libpicoquant>
- **Open Microscopy's BioFormat** supports PicoQuant bin files ([pre-histogrammed Image Data File](#)) <https://www.openmicroscopy.org/site/support/bio-formats5/formats/picoquant-bin.html>
- **dspFRET:** Matlab code for analyzing TTTR data from PicoQuant. Typical application is analysis of single pair FRET measurements with pulsed interleaved excitation (PIE-FRET). <https://github.com/RobinsonLab/dspFRET>
- **FRETbursts:** software for burst analysis of freely-diffusing single-molecule FRET (smFRET) measurements including μ s-ALEX. [FRETbursts Homepage](#) and [paper](#). Have a look also at the **PyBroMo** [simulator for single molecule FRET experiments of freely diffusing particles](#).
- **FoCuS-point:** FoCuS-point utilises advanced time-correlated single-photon counting (TCSPC) correlation algorithms along with time-gated filtering and innovative data visualisation. [FoCuS Homepage](#) and [paper](#)
- [TCSPC curve fitting for phd files \(Python\)](#)

Data Acquisition

- qudi: A modular laboratory experiment management suite which supports PicoHarp300. <https://github.com/Ulm-IQO/qudi>

Controlling PicoQuant Devices

Sepia II and Solea

- <https://github.com/PicoQuant/PQLaserDrv>

TimeHarp 260

- <https://github.com/PicoQuant/TH260-Demos>

PicoHarp 300

- PicoHarp 300 TCSPC system hardware plug-in for ScopeFoundry (Python)
- Python binding, Tango device server and Taurus gui to access the PicoHarp300 photon counter
- <https://github.com/PicoQuant/PH300-v2.x-Demos>
- <https://github.com/PicoQuant/PH300-v3.x-Demos>

HydraHarp 400

- Control HydraHarp400 (python)
- <https://github.com/PicoQuant/HH400-v1.x-Demos>
- <https://github.com/PicoQuant/HH400-v3.x-Demos>

Other

- **pyPL** wxpython integration of software and hardware for micro-photoluminescence measurements: WinSpec, PicoHarp, and Thorlabs/APT piezos and stepper motors.
<https://github.com/kaseyrussell/pyPL>
- PicoHarp Library in Python
<https://github.com/srgblnch/PicoHarp300>

GitHub repositories mentioning PicoQuant and supporting the PicoQuant t3r data file format:

<https://github.com/search?&q=picoquant+OR+picoharp+OR+hydraharp+OR+TimeHarp>

Copyright of this document belongs to PicoQuant GmbH. No parts of it may be reproduced, translated or transferred to third parties without written permission of PicoQuant GmbH. All information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuracies or omissions. Specifications and external appearances are subject to change without notice.



PicoQuant GmbH
Rudower Chaussee 29 (IGZ)
12489 Berlin
Germany

P +49-(0)30-1208820-89
F +49-(0)30-1208820-90
info@picoquant.com
www.picoquant.com