

Lectures

This page is a collection of lectures found on youtube that explain the basics of **Wfluorescence**, **fluorescence**, **fluorescen**

Fluorescence and Fluorescence Microscopy

Fluorescence is the emission of light by a substance that has absorbed light or other electromagnetic radiation. It is a form of luminescence. In most cases, the emitted light has a longer wavelength. (WFluorescence)

Fluorescence and Fluorescence Microscopy explained by Nico Stuurman at UCSF.

http://www.ibiology.org/ibioseminars/techniques/nico-stuurman-part-1.html

Confocal Microscopy

Mark Armitage explains major components of confocal microscopes and their use.

This video is part of a series of videos. This is the related playlist: https://www.youtube.com/playlist?list=PLO0yDSkagDljzI0bdM5m1njrdLdVa80 N

Nikons website MicroscopyU.com has a very nice course on confocal microscopy.

Kurt Thorn: Optical Sectioning and Confocal Microscopy

Confocal Microscopy explained by Kurt Thorn of the Nikon imaging center at UCSF

http://www.ibiology.org/ibioseminars/techniques/kurt-thorn-part-1.html

Kurt Thorn: Two Photon Microscopy

Kurt Thorn of UCSF talks about two-photon microscopy which uses intense pulsed lasers to image deep into biological samples. It can be used for imaging thick tissue specimens or even imaging inside of live animals.

http://www.ibiology.org/ibioeducation/taking-courses/two-photon-microscopy.html

Philippe Bastiaens: Fluorescence Lifetime Microscopy (FLIM)

Fluorescence-lifetime imaging microscopy or FLIM is an imaging technique for producing an image based on the

differences in the exponential decay rate of the fluorescence from a fluorescent sample. It can be used as an imaging technique in confocal microscopy, two-photon excitation microscopy, and multiphoton tomography. The lifetime of the fluorophore signal, rather than its intensity, is used to create the image in FLIM. This has the advantage of minimizing the effect of photon scattering in thick layers of sample (W FLIM).

http://www.ibiology.org/ibioeducation/taking-courses/fluorescence-lifetime-imaging-microscopy.html

Copyright of this document belongs to PicoQuant GmbH. No parts of it may be reproduced, translated or transferred to third parties without written permission of PicoQuant GmbH. All information given here is reliable to our best knowledge. However, no responsibility is assumed for possible inaccuraciesor omissions. Specifi cations and external appearances are subject to change without notice.



PicoQuant GmbH Rudower Chaussee 29 (IGZ) 12489 Berlin Germany P +49-(0)30-1208820-89 F +49-(0)30-1208820-90 info@picoquant.com www.picoquant.com